

# MINI TRAUMA

## First Aid Kit

119 Items

#FAK5036

Back pack and straps with inner detachable pouches designed for use by Doctors, Nurses, Industrial, Mining and Professional Personnel.  
For Cuts, Burns, Bleeding, IV Therapy.

Bag Dimensions: 380 mm x 160 mm x 255 mm



BURNS		
QTY	ITEM	INDICATION
4	<b>Burnshield Dressing</b> 100 mm x 100 mm	<b>Sterile trauma hydrogel</b> For the emergency treatment of first, second and third degree burns. Cools and soothes. Most effective in the first 24 hours.
1	<b>Burnshield Dressing</b> 200 mm x 200 mm	
1	<b>Burnshield Dressing</b> 600 mm x 400 mm	
2	<b>Burnshield Digit Dressing</b> 25 mm x 0,5 m	
1	<b>Burnshield Hydrogel</b> 125 ml	
BANDAGES & TAPES		
QTY	ITEM	INDICATION
2	<b>Conform Bandage</b> 75 mm x 4.5 m (Stretch)	<b>Cover and protect</b> General purpose binding. Stretch polyester.
2	<b>Conform Bandage</b> 100 mm x 4.5 m (Stretch)	
2	<b>FAD No. 2</b> 50 mm x 75 mm (Pad size)	<b>Wound pressure dressing</b> Minimises bleeding with absorbant pad.
4	<b>FAD No. 3</b> 75 mm x 100 mm (Pad size)	
4	<b>FAD No. 5</b> 150 mm x 200 mm (Pad size)	
2	<b>Trauma Pad</b> 300 mm x 200 mm	<b>Wound pressure dressing</b> Minimises bleeding on a large scale with absorbent pad.
4	<b>Triangular Bandage</b> 900 mm x 900 mm x 1270 mm	<b>Non woven</b> Arm sling which supports or immobilises an injury to a bone or joint.
2	<b>Triangular Bandage</b> Calico hemmed	<b>Re-usable</b> Arm sling which supports or immobilises an injury to a bone or joint.
4	<b>Gauze Swab</b> 75 mm x 75 mm (5's)	<b>Sterile</b> Cleaning and dressing wounds.
1	<b>Post OP Waterproof Clear</b> 50 mm x 70 mm	<b>Sterile</b> Dressing for the protection of wounds.
1	<b>Post OP Non woven + Pad</b> 90 mm x 150 mm	
1	<b>Post OP Non woven + Pad</b> 90 mm x 200 mm	
1	<b>Post OP Non woven + Pad</b> 90 mm x 60 mm	
1	<b>Fabric Plaster Roll</b> 25 mm x 3 m	<b>Adhesive fabric tape</b> Support finger sprains and strains. Securing bandages and dressings.
1	<b>Plastipore Plastic Tape</b> 25 mm x 5 m	<b>Hypo-Allergenic</b> Adhesive waterproof tape. Easy Tear. Secures dressings or bandages.
1	<b>Cohesive Soft1 Foam Bandage</b>	<b>Self Adhesive</b> Soft, highly flexible bandage. For absorbing blood and other liquids.
ACCESSORIES		
QTY	ITEM	INDICATION
4	<b>Gloves</b> (Pairs)	<b>Protective barrier</b> Prevents exposure to infection when handling open wounds.
50	<b>Plaster Strips</b>	<b>Adhesive plasters</b> Covers and protects cuts and grazes.

1	<b>Rescue Scissors</b>	Sturdy, heavy duty scissors used to cut through thick fabrics and materials.
1	<b>Emergency Rescue Blanket</b> 2 m x 1 m	<b>Thermal insulation cover</b> Reduces heat loss in the body caused by exposure. Facilitates in shock reduction.
1	<b>Eye Pad</b>	<b>Non-Sterile</b> Covers and protects injured eye.
2	<b>CPR Mouthpiece</b>	<b>Barrier device</b> Protection from exposure to infection when performing CPR.
1	<b>Airway (Guedel) 00</b>	<b>Non-sterile</b> Used to gain access to the upper airways for the administration of oxygen.
1	<b>Airway (Guedel) 2</b>	
1	<b>Airway (Guedel) 3</b>	
1	<b>Airway (Guedel) 4</b>	
1	<b>Blood Admin Set</b> 10 Drop	<b>Delivery device</b> Allows intravenous infusions of blood.
1	<b>I.V. Admin Set</b> 20 Drop	<b>Delivery device</b> Allows intravenous infusions of general fluids and medication.
1	<b>I.V. Cannula 20G</b>	<b>Delivery device</b> Needle placed into the vein for administration of blood and fluids.
1	<b>I.V. Cannula 18G</b>	
1	<b>Safety Pins (12's)</b>	Secure a bandage.
1	<b>Metal Forceps</b>	Remove splinters from a wound.
4	<b>Alcohol Pads</b>	<b>Isopropyl Alcohol 70%</b> Cleans and disinfects. Can be used on surfaces or as hand scrubs.
1	<b>Needle 21G</b>	Used for drawing blood.
1	<b>Syringe 10 ml</b>	Used for fluid withdrawal, injection of medication and wound irrigation.
1	<b>Syringe 5 ml</b>	
1	<b>Cetrimide</b> 50 ml	<b>Antiseptic wound cleaner</b> Can be applied directly on to wound.

NB: Note the items with a \* will be taken out when the kits are exported out of South Africa

**LEVTRADE INTERNATIONAL (PTY) LTD**

1 Manchester Road, Medopark Block 3, Wadeville, 1428, South Africa  
Tel: +2711440 5171 | info@levtrade.co.za | www.levtrade.co.za | Made in South Africa